

# Characterization of an i-Motif Forming Sequence from the Promoter Region of ALOX5

Ying Xia, Mark Searcey & Zoë A. E. Waller

School of Pharmacy, University of East Anglia, Norwich Research Park, Norwich, NR4 7TJ

## Abstract

i-Motifs are DNA secondary structures that form in cytosine rich sequences, consisted of four strands, stabilized by hemi-protonated cytosine–cytosine base pairs[1]. Figure 2 shows structure of i-motif. A number of i-motif forming sequences have been characterized previously[2], including measurement of their transitional pHs (pH<sub>T</sub>, the pH at which the structure is 50% folded), melting temperatures (T<sub>m</sub>, the temperature at which the structure is 50% unfolded on heating), annealing temperatures (T<sub>a</sub>, the temperature at which structure is 50% folded on cooling). Here we describe the characterization of the i-motif forming sequence from the promoter region of ALOX5, the gene encoding the enzyme arachidonate 5-lipoxygenase[3]. This gene is a current target for pharmaceutical intervention in a number of diseases, including Leukemia[4][5]. We are particularly interested in the **ALOX5** i-motif, because it has a similar sequence to other i-motifs previously characterized from the promoter regions of DAP and MSMO1 (Figure 1)[2]. This poster will describe the characterization of **ALOX5** and comparison with **DAP** and **MSMO1**. i-Motifs were assessed by circular dichroism (CD) and UV spectroscopy. We found that increasing C-tract repeats, increases thermal stability, but also complexity in the structures possible.

Figure 1 : Sequences used in this study

ALOX5	(CCCCCG) <sub>5</sub> CCCCC	n=6
DAP	(CCCCCG) <sub>4</sub> CCCCC	n=5
MSMO1	(CCCCCG) <sub>3</sub> CCCCC	n=4





Figure 2 : Structures of C-C base pairs and hTelo i-motif DNA

## Biophysical Techniques

- CD spectroscopy is used to measure transitional pH of oligonucleotide and melting temperature
- UV spectroscopy is used to measure melting temperature & annealing temperature and thermal difference spectrum

Thermal difference spectrum of ALOX5




The positive peak is at 240 nm, the negative peak is at 295 nm, It means formation of i-motif is occurring.


BLAST

ALOX5	(CCCCCG) <sub>5</sub> CCCCC	nothing	(GGGGGC) <sub>5</sub> GGGGG	nothing
DAP	(CCCCCG) <sub>4</sub> CCCCC	Setaria viridis cultivar hits	(GGGGGC) <sub>4</sub> GGGGG	Setaria viridis hits
MSMO1	(CCCCCG) <sub>3</sub> CCCCC	Setaria viridis hits	(GGGGGC) <sub>3</sub> GGGGG	Diutina rugosa hits


## CD Titrations at Different pHs




## CD Melt at pH 7.0




## 2.5 uM ALOX5 UV melt and anneal at pH7.0



## 2.5 uM ALOX5 UV melt and anneal at pH7.0 repeat



Name	T <sub>m</sub> (pH 7.0) Measured by UV	T <sub>a</sub> (pH 7.0) Measured by UV	T <sub>m</sub> (pH 7.0) Measured by CD	Transitional pH	sequence	sequence length	number of (CCCCCG) repeats
ALOX5	28.7±0.2	22.5 ± 0.2	29.97 ± 0.23	6.95	(CCCCCG) <sub>5</sub> CCC-CC	35	5
DAP	24.7 ± 0.5	22.0 ± 0.4	26.4 ± 0.36	7.1	(CCCCCG) <sub>4</sub> CCC-CC	29	4
MSMO1	16.6 ± 0.5	15.9 ± 0.4	14.7 ± 1.6	6.71	(CCCCCG) <sub>3</sub> CCC-CC	23	3



## Conclusion and future work

- Through current results of transitional pHs and T<sub>m</sub>s and T<sub>a</sub>s, we found the more C-tract repeats in sequence, the more stable in thermal stability.
- Future work will continue to complete UV absorbance experiment to get annealing temperature of ALOX5 and compare results with DAP and MSMO1.

## References:

- Gehring, K., Leroy, J.-L. & Guéron, M. *Nature* **363**, 561–565 (1993).
- Wright, E. P., Huppert, J. L. & Waller, Z. A. E. *Nucleic Acids Res.* **45**, 2951–2959 (2017).
- Funk, C. D., Hoshiko, S., Matsumoto, T., Radmark, O. & Samuelsson, B. *Proc. Natl. Acad. Sci. U. S. A.* **86**, 2587–2591 (1989).
- Ochs, M. J., Suess, B. & Steinhilber, D. *Basic Clin. Pharmacol. Toxicol.* **114**, 78–82 (2014).
- Dishart, D. et al. *Biochim. Biophys. Acta - Mol. Cell Biol. Lipids* **1738**, 37–47 (2005).